

DNA NEW PLOT SETTINGS IN ID-X

- A.1 **TOOLS>GENEMAPPER ID-X MANAGER (CTRL+M)**
- A.2 Select the **PLOT SETTINGS TAB**
- A.3 Select **NEW...** at the bottom of the dialog box
- A.4 In the **GENERAL TAB**:
 - A.4.1 Name: name table setting
 - A.4.2 Security group: Genemapper *ID-X* security group
 - A.4.3 Description: leave blank
- A.5 **SAMPLE HEADER TAB**: set to preference
- A.6 **GENOTYPES HEADER TAB**: set to preference
- A.7 **SIZING TABLE TAB**: set to preference
- A.8 **LABELS TAB**: Set to lab specifics
 - A.8.1 Label 1: allele call
 - A.8.2 Label 2: height
 - A.8.3 Label 3: data point
 - A.8.4 Labels include: assigned allele, custom allele, allelic ladder and artifacts
- A.9 **DISPLAY SETTINGS TAB**:
 - A.9.1 Remove **MARKER INDICATORS** Font Settings: set desired font
 - A.9.2 Change all other settings to meet preferences

Document ID	Revision	Approval	Date Published
1814	8	Supervising Criminalist - Biology	12/20/2018 9:33:35 AM